

ALKBH4 Rabbit mAb

ABZ1796


Product Overview

Name	ALKBH4 Rabbit mAb
Accession(Primary)	Q9NXW9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q9NXW9)

Gene Name ALKBH4

About Gene Symbol



Western blot analysis of ALKBH4 in K562, HeLa lysates using ALKBH4 antibody.

