

# ANKRD20A1 Rabbit pAb

ABZ1631

## Product Overview

<b>Name</b>	ANKRD20A1 Rabbit pAb
<b>Catalog #</b>	ABZ1631
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	Q5TYW2
<b>Application Note (Approx.)</b>	WB1:500-1:2000
<b>Precautions</b>	

## Target information(Q5TYW2)

**Synonyms****Gene ID** 84210**Other Names****Function****Cellular location****Note**



Western blot analysis of ANKRD20A1 in 3T3 lysates using ANKRD20A1 antibody.