

ANTI-AIFM1 (AIF) ANTIBODY (ABP1063)

Rabbit Polyclonal Antibody


ABP1063

Product Overview

Name	ANTI-AIFM1 (AIF) ANTIBODY (ABP1063)
Catalog #	ABP1063
Clonality	Polyclonal
Accession(Primary)	O95831
Application Note (Approx.)	WB1:1000
Precautions	

Target information(O95831)

Synonyms**Gene ID** 9131**Other Names****Function****Cellular location****Note**



Sample: Human (293) cell lysate. PVDF membrane was probed with Anti-AIFM1/AIF Antibody (Catalog# PP1063-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.