

ANTI-CLATHRIN HEAVY CHAIN ANTIBODY (ABP1024)

Rabbit Polyclonal Antibody

ABP1024

Product Overview

Name	ANTI-CLATHRIN HEAVY CHAIN ANTIBODY (ABP1024)
Catalog #	ABP1024
Clonality	Polyclonal
Accession(Primary)	Q00610
Application Note (Approx.)	WB1:1000
Precautions	

Target information(Q00610)

Synonyms


Gene ID 1213

Other Names

Function

Cellular location

Note



Sample: Human (HeLa) cell lysate. PVDF membrane was probed with Anti-Clathrin Heavy Chain Antibody (Catalog# PP1024-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.