

# Acetyl Coenzyme A Carboxylase Rabbit mAb


ABZ2197

## Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Acetyl Coenzyme A Carboxylase Rabbit mAb     |
| <b>Accession(Primary)</b>   | Q13085                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                            |
| <b>Calculated MW</b>        | Calculated MW: 266 kDa; Observed MW: 266 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:100-1:200             |


## Target information (Q13085)

|                          |                                                                                                                                                                                                                                                                                     |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | ACACA/ACACB                                                                                                                                                                                                                                                                         |
| <b>About Gene Symbol</b> | ACACA (Acetyl-CoA Carboxylase Alpha) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |




Western blot analysis of Acetyl Coenzyme A Carboxylase in A431 lysates using Acetyl Coenzyme A Carboxylase antibody.


Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of Acetyl Coenzyme A Carboxylase in C2C12 lysates using Acetyl Coenzyme A Carboxylase antibody.



Western blot analysis of AcetylCoA Carboxylase in A431 lysates using Acetyl Coenzyme A Carboxylase antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using AcetylCoA Carboxylase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.