

Acetyl-Histone H2B (Lys20) Rabbit mAb


ABZ5257

Product Overview


Name	Acetyl-Histone H2B (Lys20) Rabbit mAb
Catalog #	ABZ5257
Clonality	Monoclonal
Accession(Primary)	P3377
Application Note (Approx.)	WB1:1000-1:5000 IHCP-1:50-1:100
Precautions	

Target information(P3377)

Synonyms**Gene ID** 3018**Other Names****Function****Cellular location****Note**



Western blot analysis of Histone H2B (acetyl K20) in HeLa treated with TSA lysates using Acetyl-Histone H2B (Lys20) antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using Histone H2B (acetyl K20) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.