

# Acetyl-Histone H3 (Lys14) Rabbit mAb


ABZ4315

## Product Overview

<b>Name</b>	Acetyl-Histone H3 (Lys14) Rabbit mAb
<b>Catalog #</b>	ABZ4315
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P6843
<b>Application Note (Approx.)</b>	WB1:500-1:1000
<b>Precautions</b>	

## Target information(P6843)

**Synonyms****Gene ID** 8350**Other Names****Function****Cellular location****Note**



Immunocytochemistry analysis of Acetyl-Histone H3 (Lys14) (green) in HeLa using Acetyl-Histone H3 (Lys14) antibody, and DAPI (blue).

Western blot analysis of Acetyl-Histone H3 (Lys14) in PC-12 lysates using Acetyl-Histone H3 (Lys14) antibody. Observed band size: 15kDa.



Western blot analysis of Acetyl-Histone H3 (Lys14) in A549, HL-60 lysates using Acetyl-Histone H3 (Lys14) antibody

