

# Acetyl-p53 (Lys370) Rabbit mAb


ABZ4322

## Product Overview

<b>Name</b>	Acetyl-p53 (Lys370) Rabbit mAb
<b>Catalog #</b>	ABZ4322
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P0463
<b>Application Note (Approx.)</b>	WB1:500-1:1000
<b>Precautions</b>	

## Target information(P0463)

**Synonyms****Gene ID** 7157**Other Names****Function****Cellular location****Note**



Immunocytochemistry analysis of Acetyl-p53 (Lys370) (green) in 3T3 using Acetyl-p53 (Lys370) antibody, and DAPI(blue).

Western blot analysis of Acetyl-p53 (Lys370) in Hela, Hela+TSA, 3T3, 3T3+TSA lysates using Acetyl-p53 (Lys370) antibody.



Western blot analysis of Acetyl-p53 (Lys370) in rat Brain lysates using Acetyl-p53 (Lys370) antibody

