

## Acetyl-p53 (Lys382) (7E1) Mouse mAb


ABZ4321

### Product Overview


|                             |                                     |
|-----------------------------|-------------------------------------|
| <b>Name</b>                 | Acetyl-p53 (Lys382) (7E1) Mouse mAb |
| <b>Accession(Primary)</b>   | P04637                              |
| <b>Host Species</b>         | Mouse                               |
| <b>Clonality</b>            | Monoclonal                          |
| <b>Isotype</b>              | IgG1                                |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                   |
| <b>Recommended Dilution</b> | IHC-1:50-1:100                      |

### Target information (P04637)


|                          |                                                                                                                                                                                                                                      |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | TP53                                                                                                                                                                                                                                 |
| <b>About Gene Symbol</b> | TP53 (Tumor Protein P53) is a tumor suppressor protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl-p53 (Lys382) (7E1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Acetyl-p53 (Lys382) (7E1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using Acetyl-p53 (Lys382) (7E1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.