

Activin Receptor Type IIB Rabbit mAb

ABZ2148


Product Overview

Name	Activin Receptor Type IIB Rabbit mAb
Accession(Primary)	Q13705
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (Q13705)

Gene Name ACVR2B

About Gene Symbol



Western blot analysis of Activin Receptor Type IIB/ACVR2B in CHO-K1, mouse Brain lysates using Activin Receptor Type IIB antibody.

