


# Aldehyde dehydrogenase 10 Mouse mAb

ABZ3947

## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Aldehyde dehydrogenase 10 Mouse mAb        |
| <b>Accession(Primary)</b>   | P51648                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG2a                                      |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 55 kDa; Observed MW: 55 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

## Target information (P51648)

**Gene Name** ALDH3A2**About Gene Symbol**

Western blot analysis of Aldehyde dehydrogenase 10 in A549, U251 and HepG2 lysates using Aldehyde dehydrogenase 10 antibody.

