

## Aldolase Rabbit mAb


ABZ4983

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Aldolase Rabbit mAb                        |
| <b>Accession(Primary)</b>   | P04075                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 39 kDa; Observed MW: 39 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; FC-1:20-1:50             |

### Target information (P04075)

|                          |                                                                                                                                                                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | ALDOA                                                                                                                                                                                                                                         |
| <b>About Gene Symbol</b> | ALDOA (Aldolase, Fructose-Bisphosphate A) is a metabolic enzyme involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Aldolase in A549 lysates  
using ALDOA antibody.