

# Angiopoietin Rabbit pAb


ABZ1344

## Product Overview

<b>Name</b>	Angiopoietin Rabbit pAb
<b>Catalog #</b>	ABZ1344
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	P03950
<b>Application Note (Approx.)</b>	WB1:500-1:2000 IHC1:100-1:300
<b>Precautions</b>	

## Target information(P03950)


**Synonyms****Gene ID** 284**Other Names****Function****Cellular location****Note**



Western blot analysis of Angiopoietin 1 in A549, K562 lysates using Angiopoietin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human prostaticcancer using Angiopoietin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human prostaticcancer using Angiopoietin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.