

B7H4 Rabbit mAb


ABZ2014

Product Overview

Name	B7H4 Rabbit mAb
Accession(Primary)	Q7Z7D3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 31 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200

Target information (Q7Z7D3)

Gene Name	VTCN1
About Gene Symbol	



Western blot analysis of B7H4 in MDA-MB-468,
T47D lysates using B7H4 antibody.