


Bovine Serum Albumin Mouse mAb

ABZ4017

Product Overview

Name	Bovine Serum Albumin Mouse mAb
Accession(Primary)	P02769
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Chicken
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P02769)

Gene Name Alb**About Gene Symbol**

Western blot analysis of Bovine Serum Albumin (10G8) in Bovine Serum Albuminlysates using Bovine Serum Albumin (10G8) antibody.

