

# Bovine Serum Albumin Rabbit pAb

ABZ1488

## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Bovine Serum Albumin Rabbit pAb            |
| <b>Accession(Primary)</b>   | P02769                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Chicken                                    |
| <b>Calculated MW</b>        | Calculated MW: 69 kDa; Observed MW: 66 kDa |
| <b>Recommended Dilution</b> | WB-1:2000-1:5000                           |

## Target information (P02769)

**Gene Name** Alb

**About Gene Symbol**

