

CCR1 Rabbit pAb

ABZ1589


Product Overview

Name	CCR1 Rabbit pAb
Accession(Primary)	P32246
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (P32246)

Gene Name CCR1

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human small intestine using CCR1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using CKR 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.