


CD133 (6H10) Mouse mAb

ABZ4053


Product Overview

Name	CD133 (6H10) Mouse mAb
Accession(Primary)	O43490
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 97 kDa; Observed MW: 133 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O43490)

Gene Name PROM1**About Gene Symbol**

Western blot analysis of CD133 in CaCo2 lysates using CD133 antibody.



Western blot analysis of CD133 (6H10) in CHO-K1 cells nontransfected (A) or transfected (B) with CD133 using CD133 antibody.