


## CD16 Mouse mAb

ABZ4099


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | CD16 Mouse mAb                             |
| <b>Accession(Primary)</b>   | P08637                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 29 kDa; Observed MW: 45 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100            |

### Target information (P08637)

**Gene Name** FCGR3A**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human liver tissue using CD16 (9G6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Western blot analysis of CD16 (9G6) in Jurkat, K562 lysates using CD16 (9G6) antibody