

CD276 Rabbit mAb


ABZ1978

Product Overview

Name	CD276 Rabbit mAb
Accession(Primary)	Q8VE98
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 34 kDa; Observed MW: 90-110 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:500-1:1000; FC-1:1000-1:5000; IP-1:20-1:100

Target information (Q8VE98)

Gene Name	Cd276
About Gene Symbol	



Western blot analysis of CD276 in HepG2 lysates
using CD276 antibody.