

## CD7 Rabbit pAb


ABZ1541

### Product Overview


|                             |                                                                    |
|-----------------------------|--------------------------------------------------------------------|
| <b>Name</b>                 | CD7 Rabbit pAb                                                     |
| <b>Accession(Primary)</b>   | P09564                                                             |
| <b>Host Species</b>         | Rabbit                                                             |
| <b>Clonality</b>            | Polyclonal                                                         |
| <b>Isotype</b>              | IgG                                                                |
| <b>Species Reactivity</b>   | Human, Mouse                                                       |
| <b>Calculated MW</b>        | Calculated MW: 25 kDa; Observed MW: 25 kDa                         |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC:-1:100-1:300;<br>IF-1:50-1:200; ELISA-1:20000 |

### Target information (P09564)


|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | CD7 |
| <b>About Gene Symbol</b> |     |



Western blot analysis of CD7 in K562 lysates using CD7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human muscle using CD7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human liver using CD7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.