

ALCAM (CD166) Human

Activated Leukocyte Cell Adhesion Molecule (CD166) Human Recombinant

CDA0001

Product Overview

Name	ALCAM (CD166) Human
Accession(Primary)	Q13740
Description	Activated Leukocyte Cell Adhesion Molecule (CD166) Human Recombinant

Target information (Q13740)

Gene Name

About Gene Symbol

ALCAM (Activated Leukocyte Cell Adhesion Molecule) is a cell surface marker involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.