

BCAM Human

Basal Cell Adhesion Molecule (CD239) Human Recombinant

CDA0003

Product Overview

Name	BCAM Human
Accession(Primary)	P50895
Description	Basal Cell Adhesion Molecule (CD239) Human Recombinant

Target information (P50895)

Gene Name	
About Gene Symbol	BCAM (Basal Cell Adhesion Molecule (Lutheran Blood Group)) is a cell surface marker involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.