

CD14 Mouse, Sf9

CD14 Mouse Recombinant, Sf9

CDA0010

Product Overview

Name	CD14 Mouse, Sf9
Accession(Primary)	P10810
Immunogen	Anti-human CD14 mAb, clone PAT87H7AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CD14 protein 20-349 amino acids purified from E. coli.
Description	CD14 Mouse Recombinant, Sf9

Target information (P10810)

Gene Name

About Gene Symbol