

CD44 Human, His

CD44, His Tag Human Recombinant

CDA0080

Product Overview

Name	CD44 Human, His
Accession(Primary)	P16070
Immunogen	CD44 enriched mouse LN T cells.
Description	CD44, His Tag Human Recombinant

Target information (P16070)

Gene Name

About Gene Symbol

CD44 (CD44 Molecule (Indian Blood Group)) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify t

Background

CD44 is a cell-surface glycoprotein which take part in migration, cell-cell interactions and cell adhesion. CD44 is a receptor of hyaluronic acid which interacts with other ligands such as collagens, osteopontin and matrix metalloproteinases. CD44 has various of functions such as lymphocyte activation, tumor metastasis, recirculation & homing and hematopoiesis.