

FCGR1A Human

CD64 Human Recombinant

CDA0155

Product Overview

Name	FCGR1A Human
Accession(Primary)	P12314
Immunogen	Anti-human FCGR3A mAb, clone PFCGR3AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human FCGR3A protein 18-208 amino acids purified from E. coli.
Description	CD64 Human Recombinant

Target information (P12314)

Gene Name

About Gene Symbol