


## CKMT1 (1A6) Mouse mAb

ABZ4164

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | CKMT1 (1A6) Mouse mAb                      |
| <b>Accession(Primary)</b>   | P12532                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Human, Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 47 kDa; Observed MW: 47 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (P12532)

**Gene Name** CKMT1A**About Gene Symbol**

Western blot analysis of CKMT1 in mouse Colon, 293 and A431 lysates using CKMT1 antibody.

