

Acrp30 (108-244) Human

Adiponectin (108-244 a.a.) Human Recombinant

CYK0008

Product Overview

Name	Acrp30 (108-244) Human
Accession(Primary)	Q15848
Immunogen	Anti-human adiponectin mAb, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized with Recombinant human adiponectin amino acids 15-244 purified from E. coli.
Description	Adiponectin (108-244 a.a.) Human Recombinant

Target information (Q15848)

Gene Name**About Gene Symbol**