

AITRL Human, T7

AITRL Human Recombinant, T7 Tag

CYK0028

Product Overview

Name	AITRL Human, T7
Accession(Primary)	Q9UNG2
Immunogen	Anti-human AK1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human AK1 amino acids 1-194 purified from E. coli.
Description	AITRL Human Recombinant, T7 Tag

Target information (Q9UNG2)

Gene Name**About Gene Symbol**