

ANGPTL3 (243-460) Human

Angiopoietin-like Protein 3 (243-460 a.a.) Human Recombinant

CYK0033

Product Overview

Name	ANGPTL3 (243-460) Human
Accession(Primary)	Q9Y5C1
Immunogen	Anti-human ANGPTL3 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ANGPTL3 protein 243-460 amino acids purified from E. coli.
Description	Angiopoietin-like Protein 3 (243-460 a.a.) Human Recombinant

Target information (Q9Y5C1)

Gene Name**About Gene Symbol**