

## **APOM Human, HEK**

*Apolipoprotein-M Human Recombinant, HEK*

CYK0059

### **Product Overview**

|                           |                                                                                                                                                                                               |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>               | APOM Human, HEK                                                                                                                                                                               |
| <b>Accession(Primary)</b> | O95445                                                                                                                                                                                        |
| <b>Immunogen</b>          | Anti-human APP mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human APP amino acids 18-289 purified from E. coli. |
| <b>Description</b>        | Apolipoprotein-M Human Recombinant, HEK                                                                                                                                                       |

### **Target information (O95445)**

**Gene Name****About Gene Symbol**