

# Endoglin Human, His

*Endoglin Human Recombinant, His-Tag*

CYK0141

## Product Overview

|                           |                                     |
|---------------------------|-------------------------------------|
| <b>Name</b>               | Endoglin Human, His                 |
| <b>Accession(Primary)</b> | P17813                              |
| <b>Description</b>        | Endoglin Human Recombinant, His-Tag |

## Target information (P17813)

### Gene Name

### About Gene Symbol

ENG (Endoglin) is a cell surface marker involved in Hedgehog signaling and adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.