

Flt3 Ligand Human

Flt3 Ligand Human Recombinant

CYK0152

Product Overview

Name	Flt3 Ligand Human
Accession(Primary)	P49771
Description	Flt3 Ligand Human Recombinant

Target information (P49771)

Gene Name

About Gene Symbol

References for protein:

1. Title :Identification of Y589 and Y599 in the juxtamembrane domain of Flt3 as ligand-induced autophosphorylation sites involved in binding of Src family kinases and the protein tyrosine phosphatase SHP2. Publication :Published online before printMay 9, 2006, doi:10.1182/blood-2005-07-008896 Submitted: July 14, 2005Accepted: April 26, 2006 Blood September 1 , 2006vol. 108no. 51542-1550. Link : Flt3 Ligand prospec publication 2. Title :A role of Gab2 association in Flt3 ITD mediated Stat5 phosphorylation and cell survival. Publication :Received 2 February 2009; accepted for publication 30 March 2009 Correspondence: Lars Ronnstrand,Experimental Clinical Chemistry, WallenbergLaboratory, Entr. 46, 5th fl, Malmo" University Hospital, SE-20502 Malmo", Sweden.E-mail: lars.ronnstrand@med.lu.se;http://www.expkikemi.mas.lu.se Link : Flt3 Ligand Human prospec publication