

IL 3 Human

Interleukin-3 Human Recombinant

CYK0315

Product Overview

Name	IL 3 Human
Accession(Primary)	P08700
Description	Interleukin-3 Human Recombinant

Target information (P08700)

Gene Name

About Gene Symbol

References for protein:

1. Title :Interleukin-10 overexpression in macrophages suppresses atherosclerosis in hyperlipidemic mice. Publication :Published online before print March 30, 2010, doi: 10.1096/fj.09-148155 August 2010 The FASEB Journal vol. 24 no. 8 2869-2880 . Link : IL 3 prospec publication 2. Title :A Method for the Selection of Human Embryonic Stem Cell Sublines with High Replating Efficiency After Single-Cell Dissociation. Publication :Article first published online: 24 AUG 2006 DOI:10.1634/stemcells.2005-0657 Copyright © 2006 AlphaMed Press Link : IL 3 Human prospec publication 3. Title: Selective Activation of Human Dendritic Cells by OM-85 through a NF-kB and MAPK Dependent Pathway Publications: PloS one 8.12 (2013): e82867 Link: Interleukin-3 prospec publication 4. Title: Novel markers of normal and neoplastic human plasmacytoid dendritic cells Publication: Blood 111.7 (2008): 3778-3792. Link: Interleukin-3 Human prospec publication