

## IL 3 Mouse

*Interleukin-3 Mouse Recombinant*

CYK0318

### Product Overview

|                           |                                 |
|---------------------------|---------------------------------|
| <b>Name</b>               | IL 3 Mouse                      |
| <b>Accession(Primary)</b> | P01586                          |
| <b>Description</b>        | Interleukin-3 Mouse Recombinant |

### Target information (P01586)

#### Gene Name

#### About Gene Symbol

#### *References for protein:*

1. Title : Direct binding of Cbl to Y568 and Y936 of the stem cell factor receptor/c-Kit is required for ligand-induced ubiquitination, internalization and degradation. Publication :Biochemical Journal Immediate Publication. Published on 19 Jun 2006 as manuscript BJ20060464. Link : IL 3 Mouse prospec publication 2. Title :enhanced self -renewal of hematopoietic stem/progenitor cells mediated by the stem cell gene sall4. Publication : Link : Interleukin-3 Mouse prospec publication 3 . Title :Identification of Y589 and Y599 in the juxtamembrane domain of Flt3 as ligand-induced autophosphorylation sites involved in binding of Src family kinases and the protein tyrosine phosphatase SHP2. Publication :Published online before printMay 9, 2006, doi:10.1182/blood-2005-07-008896 Submitted: July 14, 2005Accepted: April 26, 2006 Blood September 1, 2006vol. 108no. 51542-1550. Link : <http://bloodjournal.hematologylibrary.org/content/108/5/1542.full>