

IL10 Human

Interleukin-10 Human Recombinant

CYK0357

Product Overview

Name	IL10 Human
Accession(Primary)	P22301
Description	Interleukin-10 Human Recombinant

Target information (P22301)

Gene Name

About Gene Symbol

References for protein:

1. Title :Activated signal transducer and activator of transcription-3 (STAT3) is a poor regulator of tumour necrosis factor- production by human monocytes. Publication :Article first published online: 9 JAN 2007 DOI:10.1111/j.1365-2249.2006.03291.x Link : IL10 prospec publication 2. Title :Suppressor of cytokine signalling-3 at pathological levels does not regulate lipopolysaccharide or interleukin-10 control of tumour necrosis factor-? production by human monocytes. Publication :Article first published online: 11 MAY 2006 DOI: 10.1111/j.1365-2567.2006.02383.x Link : IL10 Human prospec publication 3. Title :The Tyrosine Sulfate-rich Domains of the LRR Proteins Fibromodulin and Osteoadherin Bind Motifs of Basic Clusters in a Variety of HPR-binding Proteins, Including Bioactive Factors*. Publication :First Published on August 21, 2009, doi: 10.1074/jbc.M109.047076 October 16, 2009 The Journal of Biological Chemistry, 284, 28543-28553. Link : Interleukin-10 prospec publication