

IL12 Mouse, Sf9

Interleukin 12, Sf9 Human Recombinant

CYK0366

Product Overview

Name	IL12 Mouse, Sf9
Accession(Primary)	P43432
Immunogen	Anti-human IL12B mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human IL12B amino acids 23-328 purified from E. coli.
Description	Interleukin 12, Sf9 Human Recombinant

Target information (P43432)

Gene Name**About Gene Symbol**