

IL15 Human, HEK

Interleukin-15 Human Recombinant, HEK

CYK0372

Product Overview

Name	IL15 Human, HEK
Accession(Primary)	P40933
Immunogen	Anti-human IL15 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human IL15 protein 49-162 amino acids purified from E. coli.
Description	Interleukin-15 Human Recombinant, HEK

Target information (P40933)

Gene Name**About Gene Symbol**