

IL15RA Human

Interleukin 15 Receptor Alpha, Recombinant Human

CYK0373

Product Overview

Name	IL15RA Human
Accession(Primary)	Q13261
Immunogen	Anti-human IL16 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with Recombinant human IL16 amino acids 502-631 purified from E. coli.
Description	Interleukin 15 Receptor Alpha, Recombinant Human

Target information (Q13261)

Gene Name

About Gene Symbol