

IL31 Mouse

Interleukin-31 Mouse Recombinant

CYK0426

Product Overview

Name	IL31 Mouse
Accession(Primary)	Q6EAL8
Immunogen	Anti-human IL32 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human IL32 amino acids 1-331 purified from E. coli.
Description	Interleukin-31 Mouse Recombinant

Target information (Q6EAL8)

Gene Name**About Gene Symbol**