

IRF 3 Human

IFN Regulatory Factor-3 Human Recombinant

CYK0471

Product Overview

Name	IRF 3 Human
Accession(Primary)	Q14653
Description	IFN Regulatory Factor-3 Human Recombinant

Target information (Q14653)

Gene Name

About Gene Symbol

IRF3 (Interferon Regulatory Factor 3) is an interferon regulatory transcription factor involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.