

IRF1 Human

IFN Regulatory Factor-1 Human Recombinant

CYK0473

Product Overview

Name	IRF1 Human
Accession(Primary)	P10914
Immunogen	Anti-human IRF3 mAb, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice immunized with Recombinant human IRF3 amino acids 108-166 purified from E. coli.
Description	IFN Regulatory Factor-1 Human Recombinant

Target information (P10914)

Gene Name**About Gene Symbol**