

## TNFRSF10B Human, Sf9

*TNF Ligand Receptor Superfamily Member 10B Human Recombinant, Sf9*

CYK0578

### Product Overview

|                           |                                                                                                                                                                                                              |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>               | TNFRSF10B Human, Sf9                                                                                                                                                                                         |
| <b>Accession(Primary)</b> | O14763                                                                                                                                                                                                       |
| <b>Immunogen</b>          | Anti-human TNFRSF10C mAb, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice immunized with Recombinant human TNFRSF10C amino acids 26-236 purified from E. coli. |
| <b>Description</b>        | TNF Ligand Receptor Superfamily Member 10B Human Recombinant, Sf9                                                                                                                                            |

### Target information (O14763)

**Gene Name**

**About Gene Symbol**