

TSLP Mouse

Thymic Stromal Lymphopoietin Mouse Recombinant

CYK0617

Product Overview

Name	TSLP Mouse
Accession(Primary)	Q9JIE6
Immunogen	Anti-human TSLPR mAb, clone PAT4E7A, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TSLPR protein 23-231 amino acids purified from E. coli.
Description	Thymic Stromal Lymphopoietin Mouse Recombinant

Target information (Q9JIE6)

Gene Name**About Gene Symbol**