

TNF a Porcine

Tumor Necrosis Factor-Alpha Porcine Recombinant

CYK0664

Product Overview

Name	TNF a Porcine
Accession(Primary)	P23563
Immunogen	Anti-human TNF-a mAb, clone P10A7AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TNF-a protein 57-233 amino acids purified from E. coli.
Description	Tumor Necrosis Factor-Alpha Porcine Recombinant

Target information (P23563)

Gene Name

About Gene Symbol