

WNT3A Human, HEK

Wingless-Type MMTV Integration Site Family, Member 3A Human Recombinant, HEK

CYK0667

Product Overview

Name	WNT3A Human, HEK
Accession(Primary)	P56704
Immunogen	Anti-human YWHAB mAb is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with recombinant human YWHAB amino acids 1-246 purified from E. coli.
Description	Wingless-Type MMTV Integration Site Family, Member 3A Human Recombinant, HEK

Target information (P56704)

Gene Name

About Gene Symbol