

Cleaved-Caspase 3 p17 Rabbit pAb


ABZ4335

Product Overview


Name	Cleaved-Caspase 3 p17 Rabbit pAb
Accession(Primary)	P42574
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200

Target information (P42574)

Gene Name	CASP3
About Gene Symbol	CASP3 (Caspase 3) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cleaved-Caspase3 in various lysates using Cleaved-Caspase3 p17 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse colon tissue using Cleaved-Caspase3 p17 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.