

Collagen IV alpha 1 (2G3) Mouse mAb

ABZ4025


Product Overview

Name	Collagen IV alpha 1 (2G3) Mouse mAb
Accession(Primary)	P02462
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100


Target information (P02462)

Gene Name COL4A1

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human Liver Carcinoma Tissue using Collagen IV alpha 1 (2G3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Collagen IV alpha 1 (2G3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.