

## Cortactin Rabbit mAb


ABZ2138

### Product Overview


|                             |                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------|
| <b>Name</b>                 | Cortactin Rabbit mAb                                                               |
| <b>Accession(Primary)</b>   | Q14247                                                                             |
| <b>Host Species</b>         | Rabbit                                                                             |
| <b>Clonality</b>            | Monoclonal                                                                         |
| <b>Isotype</b>              | IgG                                                                                |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                  |
| <b>Calculated MW</b>        | Calculated MW: 62 kDa; Observed MW: 80-85 kDa                                      |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:500-1:1000; FC-1:20; IP-1:20-1:50 |

### Target information (Q14247)


|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | CTTN |
| <b>About Gene Symbol</b> |      |



Western blot analysis of Cortactin in K562, rat Brain, C6, 3T3, Hela lysates using Cortactin antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Cortactin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of  
Cortactin(green) in Hela using Cortactin  
antibody,and DAPI(blue)