

Creatine Kinase B type Rabbit mAb


ABZ5186

Product Overview

Name	Creatine Kinase B type Rabbit mAb
Accession(Primary)	P12277
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500

Target information (P12277)

Gene Name	CKB
About Gene Symbol	CKB (Creatine Kinase B) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Creatine kinase B type in SH-SY5Y lysates using Creatine Kinase B type antibody.